



Product Datasheet

CYTIP antibody (orb18666)


Description	Goat polyclonal antibody to CYTIP
Species/Host	Goat
Reactivity	Canine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Target	CYBR / PSCDBP
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 38kDa band observed in lysates of cell line MOLT4 (calculated MW of 40.0kDa according to NP_004279.3). Recommended for use at 1-3 µg/ml.
Protein Sequence	SLQRLQHSNGN
Clonality	Polyclonal
MW	40
Entrez	9595
Dilution Range	ELISA: 1:16000, WB: 0.3-1 µg/ml, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of jurkat (lane 1)...



Immunofluorescence analysis of paraforma...



Immunofluorescence analysis of paraforma...