

Product Datasheet

HMHA1 antibody (orb18670)

Description

Goat polyclonal to HMHA1

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA

Applications
Target

HMHA1 / HA-1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

To avoid confusion please note that this antibody is designed to recognise the whole protein as coded by NP_036424.2 and won't recognise the peptide found eluted from HLA-A*0201 molecules by den Haan et al (Science. 1998 Feb 13;279(5353):1054-7, PMID: 9461441)ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Preliminary experiments gave an approx 65-70kDa band in Jurkat lysate at 1 µg/ml. Please note that currently we cannot find an explanation in the literature for the band we observe given the predicted size of approx. 131kDda according to NP_036424.2. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

TLGSCRERQPEFV

Clonality

Polyclonal

MW

125

Entrez
[23526](#)
Dilution Range

ELISA: 1:32000, WB: 1ug/ml

Expiration Date

12 months from date of receipt.