


Product Datasheet

2 Cys Peroxiredoxin antibody (orb186711)

Description	2 Cys Peroxiredoxin antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human 2 Cys Peroxiredoxin (141-214/214 aa)
Target	2 Cys Peroxiredoxin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	17 kDa
Uniprot ID	P30044
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ga...