



Product Datasheet

ATF5 antibody (orb18672)

Description	Goat polyclonal antibody to ATF5
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ATF5
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 30kDa band observed in Human Heart lysates (calculated MW of 30.7kDa according to NP_036200.2). Recommended concentration: 1-3 µg/ml.
Protein Sequence	EVYKARSQRTRSC
Clonality	Polyclonal
MW	30.7
Entrez	22809
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Heart lysates.



Western blot analysis of Mouse and Rat S...