


Product Datasheet

ELF1 antibody (orb18678)


Description	Goat polyclonal antibody to ELF1
Species/Host	Goat
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ELF1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.WB: Approx 70kDa band observed in Human Tonsil, Lymph Node and in Jurkat T-cell line lysates (calculated MW of 67.5kDa according to NP_758961.1). Recommended concentration: 0.3-1.0 µg/ml.
Protein Sequence	AMKQNELLEPNSF
Clonality	Polyclonal
MW	67.5; 64.9
Entrez	1997
Dilution Range	ELISA: 1:64000, WB: 0.03-0.1 µg/ml, IHC-P: 2.5 µg/ml
Expiration Date	12 months from date of receipt.



orb18678 staining
(0.3ug/ml) of Human
To...



Western blot analysis
of Jurkat cell lys...



Immunohistochemical
staining of Human
Th...