



Product Datasheet

ARPC2 antibody (orb18679)


Description	Goat polyclonal antibody to ARPC2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Target	ARPC2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.IHC: In paraffin embedded Human Spleen shows cytoplasm staining in some splenocytes in and around a T-cell area. Recommended concentration, 3-5ug/ml.WB: Approx 33-34kDa band observed in Hela and Human Skeletal Muscle lysates (predicted MW of 34kDa according to NP_005722). Recommended for use at 0.1-0.5 µg/ml.
Protein Sequence	EMKTITGKTFSSR
Clonality	Polyclonal
MW	34.3
Entrez	10109
Dilution Range	ELISA: 1:32000, WB: 0.1-0.5 µg/ml, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



orb18679 staining (0.2ug/ml) of Hela lys...



orb18679 (3.8ug/ml) staining of paraffin...



Western blot analysis of CACO-2 (A), (0....