


Product Datasheet

ADAM17 antibody (orb18689)

Description	Goat polyclonal antibody to ADAM17
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ADAM17 / TACE
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx 100kDa band observed in human testis lysates (calculated MW of 93.0kDa according to NP_003174.3). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	LQRQNRVDSKETEC
Clonality	Polyclonal
MW	93
Entrez	6868
Dilution Range	ELISA: 1:128000, WB: 0.1-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human HeLa lysa...