


Product Datasheet

ACSL4 antibody (orb18692)


Description	Goat polyclonal antibody to ACSL4
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Target	FACL4 / ACSL4
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx 80kDa band observed in lysates of cell line HepG2 (calculated MW of 79.2kDa according to NP_075266.1). Recommended concentration: 0.1-0.3 µg/ml.
Protein Sequence	HYLKDIERYGGK
Clonality	Polyclonal
MW	74.4; 79.2
Entrez	2182
Dilution Range	ELISA: 1:128000, WB: 0.1-0.3 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



orb18692 (0.1µg/ml)
staining of HepG2
ly...



Immunohistochemical
staining of Human
Sm...



Immunohistochemical
staining of Human
Co...