


Product Datasheet

MRPL3 antibody (orb18699)

Description	Goat polyclonal antibody to MRPL3
Species/Host	Goat
Reactivity	Bovine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	MRPL3
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 35kDa band observed in lysates of cell line HepG2 (calculated MW of 38.6kDa according to NP_009139.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	NVCQPGAPSITFA
Clonality	Polyclonal
MW	38.6
Entrez	11222
Dilution Range	ELISA: 1:32000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb18699
staining
(1ug/ml) of
HepG2 cell...