

Product Datasheet

C Reactive Protein/Crp Antibody (orb18703)

Description

C Reactive Protein/Crp Antibody

Species/Host

Rabbit

Reactivity

Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Immunogen

E. coli-derived mouse Crp recombinant protein (Position: H20-S225). Mouse Crp shares 71.9% and 72.7% amino acid (aa) sequence identity with human and rat Crp, respectively.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

25360 MW

Uniprot ID


P14847

Dilution Range


Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Mouse, Rat, By Heat ELISA , 0.1-0.5µg/ml, Mouse, -Western blot, 0.1-0.5µg/ml, Rat, Mouse

Expiration Date


12 months from date of receipt.



WB analysis of Crp using anti-Crp antibody...



IHC analysis of Crp using anti-Crp antibody...



IHC analysis of Crp using anti-Crp antibody...