

Product Datasheet

SETDB2 antibody (orb18710)

Description	Goat polyclonal antibody to SETDB2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CLLD8 / SETDB2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:1000.WB: Approx 80kDa band observed in human heart lysates (calculated MW of 81.9kDa according to NP_114121.1). Recommended concentration: 0.3-1 µg/ml. This product is not suitable in WB on Mouse Heart and Rat Heart lysates.
Protein Sequence	GEKNGDAKTFWME
Clonality	Polyclonal
MW	81.9; 80.7
Entrez	83852
Dilution Range	ELISA: 1:8000, WB: 1.5-4.5 µg/ml
Expiration Date	12 months from date of receipt.



orb18710
(0.5ug/ml)
staining of
human he...