



Product Datasheet

ACTR1B antibody (orb18711)

Description	Goat polyclonal antibody to ACTR1B
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	ARP1 homolog B
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.IHC: In paraffin embedded Human Uterus shows fibrillar staining of cytoplasm in endometrial cells. Recommended concentration, 3-5ug/ml.WB: Approx 45kDa band observed in human heart and brain lysates (predicted MW of 43kDa according to NP_005726). Recommended for use at 0.5-2 µg/ml.
Protein Sequence	YEEDGSRAIHRKT
Clonality	Polyclonal
MW	42.3
Entrez	10120
Dilution Range	ELISA: 1:128000, WB: 0.5-2 µg/ml, IHC-P: 3-5µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human heart lys...



Immunohistochemical staining of Human Ut...