


Product Datasheet

PRDM1antibody (orb18713)

Description	Goat polyclonal to PRDM1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Target	PRDM1 / BLIMP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognize both reported isoforms (NP_001189.2 and NP_878911.1).ELISA: Peptide ELISA: antibody detection limit dilution 1:2000.WB: Approx 100kDa band observed in Human Brain (cerebellum) lysates (calculated MW of 91.8kDa according to NP_001189.2 and of 76.8 according to NP_878911.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	KVKQETVEPMDP
Clonality	Polyclonal
MW	91.8; 76.8
Entrez	639
Dilution Range	ELISA: 1:32000, IHC-P: 10 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Brain (ce...



Immunohistochemical staining of Human To...