



Product Datasheet

TXNRD1 antibody (orb18716)


Description	Goat polyclonal antibody to TXNRD1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	Thioredoxin Reductase 1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 60kDa band observed in Human Placenta lysates (predicted MW of 57kDa according to NP_003321). Recommended for use at 0.1-1 µg/ml.
Protein Sequence	TKRSGASILQAGC
Clonality	Polyclonal
MW	54.7; 54.7; 70.9
Entrez	7296
Dilution Range	ELISA: 1:32000, WB: 0.1-1 µg/ml, IHC-P: 2.5 µg/ml
Expiration Date	12 months from date of receipt.



orb18716 staining
(0.1ug/ml) of Human
Pl...



Immunohistochemical
staining of Human
Te...



Immunohistochemical
staining of Human
Sp...