

Product Datasheet

CHFR antibody (orb18717)

Description

Goat polyclonal antibody to CHFR

Species/Host

Goat

Reactivity

Bovine, Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC

Applications
Target

CHFR

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:64000. IHC: In paraffin embedded human brain (cerebellum) shows stained Purkinje cells. Recommended concentration, 2.5ug/ml. WB: Preliminary experiments gave an approx 80kDa band in human brain lysates at 0.5 µg/ml, this band was successfully blocked by incubation with the immunising peptide. Please note that currently we cannot find an explanation in the literature for the band we observe given the predicted size of approx. 72kDda according to NP_060693. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

KFNHICEQTRFKN

Clonality

Polyclonal

MW

73.4; 73.3; 72.0; 69.2; 63.9

Entrez


55743

Dilution Range

ELISA: 1:64000, WB: 0.5ug/ml, IHC-P: 2.5 µg/ml

Expiration Date

12 months from date of receipt.



Immunohistochemical staining of human br...