



Product Datasheet

AIRE antibody (orb18719)

Description	Goat polyclonal antibody to AIRE
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Target	AIRE
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.IHC: In paraffin embedded Human Thymic medulla shows staining of select nuclei. Recommended concentration, 2.5ug/ml.WB: Appr. 60kDa in transfected HEK293 lysates (calculated MW of 57.7kDa according to NP_000374.1). Recommended concentration: 0.01-0.03 µg/ml.Experiment Notes: ChIP: A customer has reported this antibody was successful in ChIP (no data).
Protein Sequence	QSMARPAAPFPS
Clonality	Polyclonal
MW	57.7
Entrez	326
Dilution Range	ELISA: 1:128000, WB: 1-2 µg/ml, IHC-P: 2.5 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Th...



Immunohistochemical staining of Human Th...



Western blot analysis of Human Lymph Nod...