


Product Datasheet

UBE20 antibody (orb18720)

Description	Goat polyclonal antibody to UBE20
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	E2-230K / UBE20
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Preliminary experiments in Human Colon lysates and in lysates of cell line HeLa gave no specific signal but low background (at antibody concentration up to 1 µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	ALLEAGMPECTEDK
Clonality	Polyclonal
MW	141
Entrez	63893
Dilution Range	ELISA: 1:32000, WB: 1 µg/ml, IHC-P: 4-6µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human HeLa cells