

Product Datasheet

CDCP1 antibody (orb18721)

Description

Goat polyclonal antibody to CDCP1

Species/Host

Goat

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Target

CDCP1 (isoform 1: C term)

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:32000. WB: Approx 90kDa band observed in Human Colon lysates (calculated MW of 92.3kDa according to NP_073753.3). Recommended concentration: 1-3 µg/ml. An additional band of unknown identity was also consistently observed at 200kDa. This band was successfully blocked by incubation with the immunising peptide. We would appreciate any feedback from people in the field - have any such results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

Protein Sequence

PLLNTQEPMEPAE

Clonality

Polyclonal

MW

92.3

Entrez


[64866](#)

Dilution Range

ELISA: 1:32000, WB: 1-3 µg/ml, IHC-P: 5 µg/ml

Expiration Date


12 months from date of receipt.



Western blot analysis of HEK293 cells us...



Western blot analysis of Human Colon lys...



Immunohistochemical staining of Human Pr...