

Product Datasheet

hsa-miR-3622b-5p miRNA Inhibitor (orb1872392)

Description	Synthetic miRNA inhibitors are used to inhibit hsa-miR-3622b-5p by transfection.
Reactivity	Human
Form/Appearance	Lyophilized powder
Storage	Shipped at 4 °C. Store at -20 °C for one year. Avoid freeze-thaw cycles after reconstitution.
Note	For research use only
Application notes	We recommend re-suspending the lyophilized synthetic miRNA using DNase and RNase-free ddH2O. To make a 100 uM stock, dissolve the lyophilized powder using 100 ul of ddH2O.
Source	Synthetic
Expiration Date	12 months from date of receipt.