


Product Datasheet

Syntenin 2 antibody (orb18724)

Description	Goat polyclonal to Syntenin 2
Species/Host	Goat
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	Syntenin 2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	This antibody is expected to recognise both human isoforms (NP_056500; NP_536737).ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Preliminary experiments gave no signal but low background in Human Brain and Liver lysates at up to 1 µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	LLHHTMDHSIPDA
Clonality	Polyclonal
MW	22.7, 31.6
Entrez	27111
Dilution Range	ELISA: 1:16000, WB: 1 µg/ml, IHC-P: 4-6µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Liver tissue showing Syntenin 2 expression.