


Product Datasheet

DCLRE1A antibody (orb18725)

Description	Goat polyclonal antibody to DCLRE1A
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	SNM1 (C Term)
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Preliminary experiments gave no signal but low background in Human Brain and Testis lysates at up to 1 µg/ml (note whole-cell RIPA lysates not nuclear lysates. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	EKYFREWKLEAGY
Clonality	Polyclonal
MW	116
Entrez	9937
Dilution Range	ELISA: 1:16000, WB: 1 µg/ml, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



orb18725
(5µg/ml)
staining of
paraffin ...