

Product Datasheet

DCL-1 antibody (FITC) (orb187282)

Description	Rabbit polyclonal antibody to Anti-DCL-1/FITC
Species/Host	Rabbit
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human DCL-1
Target	CD302
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20 for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	26kDa
Uniprot ID	Q8IX05
Dilution Range	IF=1:100-500,ICC=1:100-500
Expiration Date	12 months from date of receipt.