



Product Datasheet

MID2 antibody (orb18730)

Description	Goat polyclonal antibody to MID2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	MID2 / TRIM1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.IHC: In paraffin embedded Human Kidney shows strong staining in cytoplasm and brush border of some epithelial cells of proximal convoluted tubules. Recommended concentration, 5-10ug/ml.WB: Approx 80kDa band observed in mouse heart lysates (calculated MW of 81.3kDa according to NP_036348.1 and 77.9kDa according to NP_438112.1). Recommended concentration: 1-3 µg/ml.
Protein Sequence	QESPYVSGMKTCH
Clonality	Polyclonal
MW	83.2, 79.9
Entrez	11043
Dilution Range	ELISA: 1:32000, WB: 0.1-0.3 µg/ml, IHC-P: 5-10µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse and Rat H...



Immunohistochemical staining of Human Ki...