


Product Datasheet

PLEKHA8 antibody (orb18735)

Description	Goat polyclonal antibody to PLEKHA8
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	FAPP2 / PLEKHA8
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:64000.IHC: In paraffin embedded Human Placenta shows strong staining of cytoplasm in the margins of villi. Recommended concentration, 2-4ug/ml.WB: Preliminary experiments gave no signal but low background in Human Kidney and Testis lysates at up to 1 µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	DIQTALRNPTENT
Clonality	Polyclonal
MW	49.5
Entrez	84725
Dilution Range	ELISA: 1:64000, WB: 1 µg/ml, IHC-P: 2-4µg/ml
Expiration Date	12 months from date of receipt.



orb18735
(2.5ug/ml)
staining of
paraffin...