

Product Datasheet

rno-miR-3102 miRNA Mimic (orb1873570)

Description	Synthetic miRNA mimics are used to overexpress rno-miR-3102 by transfection.
Reactivity	Rat
Form/Appearance	Lyophilized powder
Storage	Shipped at 4 °C. Store at -20 °C for one year. Avoid freeze-thaw cycles after reconstitution.
Note	For research use only
Application notes	We recommend re-suspending the lyophilized synthetic miRNA using DNase and RNase-free ddH ₂ O. To make a 100 uM stock, dissolve the lyophilized powder using 50 ul of ddH ₂ O.
Source	Synthetic
Expiration Date	12 months from date of receipt.