



Product Datasheet

RPS6KB1 antibody (orb18737)


Description	Goat polyclonal antibody to RPS6KB1
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	p70S6K / RPS6KB1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx 65-70kDa band observed in Human Skeletal Muscle and Placenta lysates (predicted MW of 60kDa according to NP_003152.1). Recommended for use at 0.2-1 µg/ml.
Protein Sequence	MISKRPEHLRMNL
Clonality	Polyclonal
MW	59.1
Entrez	6198
Dilution Range	ELISA: 1:64000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Placenta ...



Immunohistochemical staining of Human Ki...



Western blot analysis of NIH3T3 cell lys...