


Product Datasheet

FAF2 antibody (orb18739)

Description	Goat polyclonal antibody to FAF2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	ETEA / UBXD8
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 55kDa band observed in lysates of Human and Mouse Thymus and of cell line Jurkat (calculated MW of 52.6kDa according to NP_055428.1). Recommended concentration: 0.3-1µg/ml.
Protein Sequence	HTEVLVQDLTDE
Clonality	Polyclonal
MW	52.6
Entrez	23197
Dilution Range	ELISA: 1:16000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



orb18739
staining
(0.3µg/ml) of
Jurkat ...



orb18739
(1µg/ml)
staining of
Human (A)...