



Product Datasheet

CORO1A antibody (orb18743)


Description	Goat polyclonal antibody to CORO1A
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC, WB
Target	Coronin 1 / TACO
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 60kDa band observed in MOLT4 and Jurkat lysates (calculated MW of 51.0 kDa according to NP_009005.1). In transfected HEK293 transiently expressing CORO1A band of approx. 50kDa is observed. This band is not observed in the non-transfected HEK293. Recommended for use at 0.3-1 µg/ml.
Protein Sequence	KRLDRLEETVQAK
Clonality	Polyclonal
MW	51
Entrez	11151
Dilution Range	ELISA: 1:16000, WB: 0.3-1 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat lysate u...



Western blot analysis of HEK293 lysate u...



Western blot analysis of Human Lymphobla...