



Product Datasheet

FYN antibody (orb18747)

Description	Goat polyclonal antibody to FYN
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	FYN
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 55-60kDa band observed in Jurkat, A431 and NIH-3T3 lysates (predicted MW of 62kDa according to NP_002028.1). Recommended for use at 0.2-1 µg/ml.
Protein Sequence	VQCKDKEATKLTE
Clonality	Polyclonal
MW	60.8; 60.1; 54.5
Entrez	2534
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat lysate U...



Western blot analysis of A431 (Lane 1) a...