


Product Datasheet

CYB561D2 antibody (orb18749)

Description	Goat polyclonal to CYB561D2
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CYB561D2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 25kDa band observed in human liver lysates (predicted MW of 26kDa according to NP_008953). Recommended for use at 1-3 µg/ml.
Protein Sequence	VSNAYLYRKRIQP
Clonality	Polyclonal
MW	24
Entrez	11068
Dilution Range	ELISA: 1:16000, WB: 1-3 µg/ml
Expiration Date	12 months from date of receipt.



orb18749
staining
(1ug/ml) of
human live...