



Product Datasheet

beta Defensin 3 antibody (FITC) (orb187581)


Description	beta Defensin 3 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Tested Applications	FC, ICC, IF
Immunogen	KLH conjugated synthetic peptide derived from human beta Defensin 3
Target	beta Defensin 3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P81534
Dilution Range	Flow-Cyt=1µg/Test,IF=1:200-800,ICC=1:100-500
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of Human kidney tissue.



Flow cytometric analysis of HepG2 Cells.



Immunofluorescence analysis of Human colon tissue.