

## Product Datasheet

### CDC42BPA antibody (orb18762)

**Description**

Goat polyclonal antibody to CDC42BPA

**Species/Host**

Goat

**Reactivity**

Bovine, Canine, Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested**

ELISA, IHC

**Applications**
**Target**

Cdc42-binding kinase alpha

**Preservatives**

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

**Storage**

Aliquot and store at -20°C. Minimize freezing and thawing.

**Note**

For research use only

**Application notes**

ELISA: Peptide ELISA: antibody detection limit dilution 1:16000. WB: Preliminary experiments gave an approx 30kDa and 40kDa bands in Human Pancreas lysates at 1 µg/ml. Please note that currently we cannot find an explanation in the literature for the band we observe given the calculated sizes of 186kDa according to NP\_055641.3 and 196kDa according to NP\_003598.2. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported

**Protein Sequence**

SGEVRLRQLEQFI

**Clonality**

Polyclonal

**MW**

186, 196

**Entrez**


8476

**Dilution Range**

ELISA: 1:16000, IHC-P: 3.75 µg/ml

**Expiration Date**

12 months from date of receipt.



Immunohistochemical staining of Human Ce...