



Product Datasheet

SMUG1 antibody (orb18767)

Description	Goat polyclonal antibody to SMUG1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Target	SMUG1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Approx 35kDa band observed in Molt-4 lysates (predicted MW of 31kDa according to NP_055126. Recommended for use at 0.5-2 µg/ml.
Protein Sequence	PQAFLLGSIHEPA
Clonality	Polyclonal
MW	29.9; 19.6; 27.0;; 17.6
Entrez	23583
Dilution Range	ELISA: 1:128000, WB: 0.1 µg/ml, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Molt-4 lysate u...



Immunohistochemical staining of Human Co...