


Product Datasheet

AKAP10 antibody (orb18769)

Description	Goat polyclonal antibody to AKAP10
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	AKAP10 / D-AKAP2
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Preliminary experiments gave no signal but low background in Human Skeletal Muscle lysates at up to 1 µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	AQYDQPLEKSTKL
Clonality	Polyclonal
MW	73.8
Entrez	11216
Dilution Range	ELISA: 1:16000, WB: 0.3-1.0ug/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of Human Testis lysate