



Product Datasheet

GRIPAP1 antibody (orb18773)

Description	Goat polyclonal antibody to GRIPAP1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Target	GRASP1 / GRIPAP1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:32000.WB: Preliminary experiments gave no signal but low background in Human Brain lysates at up to 1 µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates
Protein Sequence	LREQTTGLAAELQ
Clonality	Polyclonal
MW	96.0; 71.9
Entrez	56850
Dilution Range	ELISA: 1:32000, WB: 1 µg/ml, IHC-P: 3.75 µg/ml
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human Tissue



Immunohistochemical staining of Human Cortex