

Product Datasheet

ARHGAP26 antibody (orb18780)

Description

Goat polyclonal antibody to ARHGAP26

Species/Host

Goat

Reactivity

Bovine, Canine, Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

orb18780 (0.3ug/ml)
staining of Human
He...

Target

GRAF / ARHGAP26

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage


Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: Peptide ELISA: antibody detection limit dilution 1:64000. WB: Approx 90kDa band observed in Human Heart lysates (calculated MW of 92.2kDa according to NP_055886.1). Recommended concentration: 0.3-1 µg/ml. An additional strong band of unknown identity was also consistently observed at 16kDa. This band was successfully blocked by incubation with the immunising peptide. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates Have any further splice variants/modified forms been reported


Immunohistochemical
staining of Human
Co...

Protein Sequence

KTGLIPENYVEFL

Clonality

Polyclonal

MW

92.2; 86.2

Entrez

23092

Dilution Range

ELISA: 1:64000, WB: 0.3-1 µg/ml, IHC-P: 3.75 µg/ml

Expiration Date

12 months from date of receipt.