

Product Datasheet

hsa-miR-499a-3p miRNA Mimic (orb1878148)

Description	Synthetic miRNA mimics are used to overexpress hsa-miR-499a-3p by transfection.
Reactivity	Human
Form/Appearance	Lyophilized powder
Storage	Shipped at 4 °C. Store at -20 °C for one year. Avoid freeze-thaw cycles after reconstitution.
Note	For research use only
Application notes	We recommend re-suspending the lyophilized synthetic miRNA using DNase and RNase-free ddH ₂ O. To make a 100 uM stock, dissolve the lyophilized powder using 50 ul of ddH ₂ O.
Source	Synthetic
Expiration Date	12 months from date of receipt.