



Product Datasheet

Goat Anti-Mouse/Rabbit IgG (H&L)-HRP Polymer (Ready to use) (orb1878754)

Description	Ready to use Goat Polyclonal Secondary Antibody to Mouse/Rabbit IgG (H&L) HRP polymer labeled
Reactivity	Mouse, Rabbit
Conjugation	HRP
Tested Applications	ELISA, IH, WB
Immunogen	Mouse/Rabbit IgG
Form/Appearance	Liquid in PBS, pH 7.2, containing 1% rAlbumin, 0.05% Proclin
Storage	Shipped and store at 4°C for one year.
Note	For research use only
Clonality	Polyclonal
Source	Goat
Dilution Range	Ready to use
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis staining in...



Immunohistochemical analysis staining in...