


Product Datasheet

Goat Anti-Rabbit IgG (H&L)-HRP Polymer (orb1878758)

Description	Goat Polyclonal Secondary Antibody to Rabbit IgG (H&L) HRP polymer labeled
Reactivity	Rabbit
Conjugation	HRP
Tested Applications	ELISA, IH, WB
Immunogen	Rabbit IgG
Form/Appearance	0.5 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.05% Proclin
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Goat
Dilution Range	E (1/5000 - 1/20000), WB (1/5000 - 1/20000), IH (1/100 - 1/500), IC (1/100 - 1/500)
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis staining in...