

## Product Datasheet

### Mouse Anti-Human IgA - AP (orb1878887)

|                            |                                                                                                                      |
|----------------------------|----------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse Monoclonal Secondary Antibody to Human IgA AP labeled                                                          |
| <b>Reactivity</b>          | Human                                                                                                                |
| <b>Conjugation</b>         | AP                                                                                                                   |
| <b>Tested Applications</b> | ELISA, IH, WB                                                                                                        |
| <b>Immunogen</b>           | Human IgA                                                                                                            |
| <b>Form/Appearance</b>     | 1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                     |
| <b>Note</b>                | For research use only                                                                                                |
| <b>Clonality</b>           | Monoclonal                                                                                                           |
| <b>Source</b>              | Mouse                                                                                                                |
| <b>Dilution Range</b>      | E (1/10000 - 1/50000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                      |