

Product Datasheet

Rabbit Anti-Donkey IgG (H&L) - AP (orb1878904)

Description	Rabbit Polyclonal Secondary Antibody to Donkey IgG (H&L) AP labeled
Reactivity	Donkey
Conjugation	AP
Tested Applications	ELISA, IH, WB
Immunogen	Donkey IgG
Form/Appearance	1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Dilution Range	E (1/10000 - 1/50000), WB (1/5000 - 1/20000), IH (1/500 - 1/1000), IC (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.