



Product Datasheet

Anti-ak2 Antibody (orb1878966)


Description	Rabbit polyclonal antibody to ak2
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	FC, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-terminal region of zebrafish ak2. The exact sequence is proprietary.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q1L8L9
Entrez	321793
Dilution Range	WB (1/500 - 1/1000), IH (1/50 - 1/200), FC (1/10 - 1/30)
Expiration Date	12 months from date of receipt.



Western blot analysis of ak2 expression ...



Immunohistochemical analysis of ak2 stai...



Flow cytometric analysis of ZF4 cells us...