


Product Datasheet

Anti-USP17L24 Antibody (orb1879006)

Description	Rabbit polyclonal antibody to USP17L24
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-terminal region of human USP17L24. The exact sequence is proprietary.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q0WX57
Entrez	728369 , 728405 , 728379 , 728400 , 728393 , 728373 , 728419
Dilution Range	WB (1/500 - 1/1000), FC (1/10 - 1/50)
Expiration Date	12 months from date of receipt.



Western blot analysis of USP17L24 expres...



Flow cytometric analysis of HepG2 cells
...