

Product Datasheet

MORF4L1 antibody (orb18792)

Description

Goat polyclonal antibody to MORF4L1

Species/Host

Goat

Reactivity

Bovine, Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC, WB

Applications

Target

MRG15 / MORF4L1

Preservatives

Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Note

For research use only

Application notes

ELISA: antibody detection limit dilution 1:128000. IHC: In paraffin embedded Human Uterus shows nuclear staining in the epithelial cells of the endometrial glands. Data obtained from the previous goat. Recommended concentration, 2-4 µg/ml. WB: Approx 37kDa band observed in nuclear lysates of cell lines HepG2 and NIH3T3 and approx 40kDa band observed in Mouse, Rat and Pig Liver lysates (calculated MW of 37.2kDa according to Human NP_006782.1, Mouse NP_077751.1, Rat NP_001011999.1 and Pig XP_001927776.1 and 41.5kDa according to Human NP_996670.1, Mouse NP_001034236.1 and Pig XP_001927789.1). Recommended concentration: 0.1-0.3 µg/ml.

Protein Sequence

APKQDPKPKFQEG

Clonality

Polyclonal

MW

37.2; 41.5

Entrez


10933

Dilution Range


ELISA: 1:128000, WB: 0.1-0.3 µg/ml

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of Human Ut...



Immunohistochemical staining of Human Pl...



Western blot analysis of NIH3T3 (Lane 1)...