



## Product Datasheet

### Anti-SPCS3 Antibody (orb1879570)


|                            |                                                                                                                                          |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to SPCS3                                                                                                      |
| <b>Reactivity</b>          | Human, Mouse                                                                                                                             |
| <b>Conjugation</b>         | Unconjugated                                                                                                                             |
| <b>Tested Applications</b> | FC, IH, WB                                                                                                                               |
| <b>Immunogen</b>           | KLH-conjugated synthetic peptide encompassing a sequence within the C-terminal region of human SPCS3. The exact sequence is proprietary. |
| <b>Form/Appearance</b>     | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                                |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                                         |
| <b>Note</b>                | For research use only                                                                                                                    |
| <b>Clonality</b>           | Polyclonal                                                                                                                               |
| <b>Source</b>              | Rabbit                                                                                                                                   |
| <b>Uniprot ID</b>          | <b>P61009, Q9D365</b>                                                                                                                    |
| <b>Entrez</b>              | <b>76976, 76219, 60559</b>                                                                                                               |
| <b>Dilution Range</b>      | WB (1/500 - 1/1000), IH (1/50 - 1/200), FC (1/10 - 1/50)                                                                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                          |



Western blot analysis of SPCS3 expression...



Immunohistochemical analysis of SPCS3 staining...



Flow cytometric analysis of HepG2 cells ...