



## Product Datasheet

### Anti-KIAA0226L Antibody (orb1879601)

|                            |                                                                                                                                           |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to KIAA0226L                                                                                                   |
| <b>Reactivity</b>          | Human                                                                                                                                     |
| <b>Conjugation</b>         | Unconjugated                                                                                                                              |
| <b>Tested Applications</b> | FC, WB                                                                                                                                    |
| <b>Immunogen</b>           | KLH-conjugated synthetic peptide encompassing a sequence within the Central region of human KIAA0226L. The exact sequence is proprietary. |
| <b>Form/Appearance</b>     | Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.                                 |
| <b>Storage</b>             | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.                                          |
| <b>Note</b>                | For research use only                                                                                                                     |
| <b>Clonality</b>           | Polyclonal                                                                                                                                |
| <b>Source</b>              | Rabbit                                                                                                                                    |
| <b>Uniprot ID</b>          | <b>Q9H714</b>                                                                                                                             |
| <b>Entrez</b>              | <b>80183</b>                                                                                                                              |
| <b>Dilution Range</b>      | WB (1/500 - 1/1000), FC (1/10 - 1/50)                                                                                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                           |



Western blot analysis of KIAA0226L expression...



Flow cytometric analysis of Ramos cells  
...