


Product Datasheet

Anti-MPP7 Antibody (orb1879777)


Description	Rabbit polyclonal antibody to MPP7
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IH, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-terminal region of human MPP7. The exact sequence is proprietary.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q5T2T1
Entrez	143098
Dilution Range	WB (1/500 - 1/1000), IH (1/50 - 1/200), FC (1/10 - 1/30)
Expiration Date	12 months from date of receipt.



Western blot analysis of MPP7 expression...



Immunohistochemical analysis of MPP7 staining...



Flow cytometric analysis of HeLa cells expressing MPP7...