


Product Datasheet

Anti-PRSS57 Antibody (orb1879881)


Description	Rabbit polyclonal antibody to PRSS57
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-terminal region of human PRSS57. The exact sequence is proprietary.
Form/Appearance	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Clonality	Polyclonal
Source	Rabbit
Uniprot ID	Q6UWY2
Entrez	400668
Dilution Range	WB (1/1000 - 1/2000), IF/IC (1/10 - 1/50), FC (1/10 - 1/30)
Expiration Date	12 months from date of receipt.



Western blot analysis of PRSS57 expression...



Immunofluorescent analysis of Anti-PRSS57...



Flow cytometric analysis of A549 cells u...