


Product Datasheet

CPEB1 antibody (orb18801)

Description	Goat polyclonal antibody to CPEB1
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Target	CPEB1
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:128000.WB: Approx 65-70kDa band observed in human brain lysates (predicted MW of 65kDa according to NP_085097). Recommended for use at 0.5-2 µg/ml. Please note this band was always accompanied by a slightly stonger band at approx 48kDa and further both bands could be blocked by incubation with the immunising peptide.
Protein Sequence	PLMRNQKNRDSS
Clonality	Polyclonal
MW	62.1; 54.2; 53.6
Entrez	64506
Dilution Range	ELISA: 1:128000, WB: 0.5-2 µg/ml
Expiration Date	12 months from date of receipt.



orb18801
staining
(0.5ug/ml) of
human br...