


Product Datasheet

RAN antibody (orb18802)


Description	Goat polyclonal antibody to RAN
Species/Host	Goat
Reactivity	Bovine, Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Target	RAN
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	ELISA: Peptide ELISA: antibody detection limit dilution 1:16000.WB: Approx 25kDa band observed in Jurkat and A431 lysates (predicted MW of 25kDa according to NP_006316). Recommended for use at 1-3 µg/ml.
Protein Sequence	AAQGEPQVQFKLV
Clonality	Polyclonal
MW	24.4
Entrez	5901
Dilution Range	ELISA: 1:32000, WB: 0.3-1 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of A431 lysate using RAN antibody (orb18802).



Immunohistochemical staining of Human Toxinogen (Tox) tissue using RAN antibody (orb18802).



Immunohistochemical staining of Human Toxinogen (Tox) tissue using RAN antibody (orb18802).